

AUTOCLAVE

SAFE & EASY OPERATIONS
Autoclave for modern microbiologist




**UNIVERSAL
LAB SOLUTIONS**
"The Complete Laboratory Solution"

Model	Working Chamber Sizes
# 0421 – SLD	Ø 30 x 50 cm
# 0431 – SLD	Ø 35 x 55 cm
# 0441 – SLD	Ø 40 x 60 cm
# 0440 – SLD	Ø 45 x 60 cm
# 0451 – SLD	Ø 51 x 76 cm

- ~ Sole Locking Device [SLD] very easy to open and close for smart staff.(Especially lady microbiologist)
- ~ Low water level cut off – A water level monitoring system maintains constant water level and ensures safe operation of the heaters.
- ~ Self- regulating safety for over temperature
- ~ Easy-going drain from the chamber
- ~ Automatic purging and exhaust for SLD models
- ~ Safe to handle external body eliminating heat burns
- ~ Lid Lock Safety – A safety device prevents the operator from opening the door when the chamber is pressurized.
- ~ Safety Valves – A chamber is equipped with safety valves – if, the pressure exceeds the allowed safety limits the safety valves will discharge.

* Supplied with Test Certificates & calibration reports
along with traceability for Temperature Controller and Pressure Gauge.

** Images are just for illustration purpose. Specs are subject to change without prior notification



*Documentation Scanner Available as Option



~ 9320645189