

Table Top Autoclave Sterilizer Class B

This product is a mini table type pulse vacuum sterilizer with rapid cooling system. It is safe and auto-controlled by computer. Its sterilizing temperature range is 105 to 134°C, and timer range is 1 to 99 minutes.



Features:

4. The whole process of water adding, temperature rising, sterilizing, drying steam discharging controlled automatically.
5. Steam-water inner circulation system: no steam discharge, and the environment for sterilizing will be clean & dry.
6. Exhaust the cool air automatically.
7. Door safety lock system.
8. With four stainless steel sterilizing plates.
9. The chamber of the sterilizer is made of stainless steel.
10. Automatically shut off with beep reminding after sterilization.
11. With drying function.

Specification:

MODEL	TOP-250B20	TOP-250B25
Sterilizing chamber volume:	20L (Φ 250*420mm)	25L (Φ 250*520mm)
Maximum working pressure	0.22 Mpa	
Maximum working temperature	134°C	
Adjustment of temperature	105°C ~ 134°C	
Timer	1-99 min	
Chamber temperature Accuracy	± 1°C	
Source power	1.5KW/AC 220V 50Hz	
Sterilizing plate	340*200*30 mm (3 piece)	400*200*30 mm (3 piece)
Dimension	480*480*380mm	580*480*384mm
Package dimension	720*580*500mm	800*580*500mm
N/G. Weight	40/43 Kg	45/50 Kg