



## HOT PLATE (ROUND / RECTANGULAR)

- Body Made Of Thick M.s./s.s.
- Top Plate Over Hangs By Supporting Of Ceramic Beeds Which Ensures Air Insulation Between Plate And Base Unit.
- Selection Cast Iron / S.s.plate
- Maximum Surface Temperature 300°C.
- Temp. Controlled By Energy Regulator.
- Long Last Nichrome Elements.
- Various Sizes Available
- To Work On 230 Volts A. / C. Only.

Sizes available: (Round) Rating watts - 1500  
20 Cm (8") dia. Rectangular Rating watts

S.No.	SIZE AVAILABLE	RATING WATTS.
1.	25 x 30 cm.	150
2.	25 x 40 cm.	200
3.	30 x 35 cm.	300
4.	45 x 60 cm.	450



*Universal Lab Solutions*

402, Nandanvan Tower, Plot No. 125A/7,

Panvel - 410206, Dist.- Raigad

E-mail : [uls.instrument@gmail.com](mailto:uls.instrument@gmail.com)

Mob: +91- 9320645189